

MEDOP 102 Clinical



Respiratory protection: Filter Masks

Description and Composition:

Single use mask for hygiene purposes. Includes customizable nasal clip.

Components:

- Polyester
- Rubber Band
- Nasal clip

Maximum level of use: <1xTLV

CE



Elastic strap.



Smooth covering.



Adaptable clip.

Ref.	Product
902.342	MEDOP 102 Clinical

Characteristics		
Active Carbon	X	
Dolomite Test	X	
With Valve	X	
Latex free	X	
Metal free	X	
Individual bags	X	(
Size	Unic	
Weight	4 gr.	

Respiratory protection: Filter Masks

Standard and Certificacion	CE
Aplications	For non-toxic particles. May be used with concentrations of up to 1 TLV. Never use to filter gases or vapours. Do not use in environments containing less than 17% oxygen by volume, nor in explosive atmospheres. USES: sanitary. Prevents particle spreading and the spread of disease. Also suitable to prevent contagion in epidemics, cleaning of buildings or for allergic people at risk.
Maintenance Storage - Expiry date	Store the equipment in its package, in a cool, well-ventilated place, avoiding moisture, filth and dust.
Instructions Use - How to Employ	Dispose of the mask after use. Designed for one-time use only. Dispose of the mask if it offers great resistance to breathing.
Presentation	Box of 50 units. Carton of 200 boxes.
Barcode	GTIN-13: 8423173110833 GTIN-14: 88423173110839

