



# EYEWASH BOTTLES



· Pocket eyewash



· Emergency shower



· Wall-mounted eyewash



· Transportation eyewash

## Eye wear accessories

### Description and composition:

**0.9% sterile saline solution**, to be used when a product or substance reaches eyes.

**Polypropylene bottle**, an inert material that allows for optimal conservation of the solution.

**900438** - Refillable eyewash bottle: Emergency shower. For use as an emergency eye bath.  
Drainage valve. Water-tight seal.

**900303** - Pocket eyewash bottle: Contains 200 ml of 0.9% sterile saline solution.

**914437** - Wall-mounted eyewash bottle: With wall support. One wall support and 2 one litre bottles of 0.9% sterile saline solution.

**911797** - Transportation eyewash bottle: Special support for TRANSPORTATION and 2 one litre bottles of 0.9% sterile saline solution.

EN 15154-4:2009 EN 556-1:2001  
ISO 17665-1:2006 ISO 14644-1:2015 ISO 17141:2020 ISO 15223-1:2021

Ref.	Product	Capacity
900438	Refillable eyewash bottle: Emergency shower	775 cm <sup>3</sup>
900303	Pocket eyewash bottle	200 ml
914437	Wall-mounted eyewash bottle	2 one litre bottles
911797	Transportation eyewash bottle	2 one litre bottles
900355	Replacement eyewash bottle.	2 one litre bottles

## Eye wear accessories

<b>Standard and certification</b>	EN 15154-4:2009 EN 556-1:2001 ISO 17665-1:2006 ISO 14644-1:2015 ISO 17141:2020 ISO 15223-1:2021
<b>Applications</b>	0.9% sterile saline solution, to be used when a product or substance reaches eyes.
<b>Conservation Storage - Expiry</b>	Bottles should be stored in specially designed supports. Optimum storage temperature is 2 °C to 35 °C. Expiry date: see package.
<b>Directions Use</b>	Apply directly to eyes. Since the liquid is sterile, it will not harm the eyes. The product is easy to open. Precise product flow control. Once the bottle has been opened, the contents must be used within several hours. Any remaining liquid must be discarded.
<b>Presentation</b>	<b>Refillable eyewash bottle: Emergency shower:</b> Individually packed: 1 unit. 6 units per carton.  <b>Pocket eyewash bottles:</b> Individually packed: 1 unit. 12 units per carton.  <b>Wall-mounted eyewash bottles:</b> Individually packed: 1 unit. 12 units per carton.  <b>Transportation eyewash bottle:</b> Individually packed: 1 unit. 12 units per carton.  <b>Replacement eyewash bottle:</b> Individually packed: 1 unit. 12 units per carton.

