




# Lumbar back support belt



· Lumbar back support belt with suspenders



· Strapless lumbar back support belt

## Body protection: ergonomic line

### Description and composition:

**Lumbar back support in polyester**, with four stays to provide proper stability.

Back support made of PP (polypropylene) with Velcro straps (in nylon) for better protection of the lumbar area. The lumbar back support with straps has POM (polyoxymethylene) buckles on each strap.

Provides support for the lower back and abdomen, allowing freedom of movement and maximum comfort due to its anatomical shape.

**Has a Velcro fastener and additional adjustment straps also** (Velcro in Nylon). Meaning better fit and more mobility. This model is available with suspenders for greater comfort and to avoid losing the belt.

**Washable and breathable.** Dynamic, flexible fastening which does not prevent movement. No mechanical parts that can break or need changing. Includes four stays.


**Size:** S, M, L, XL  
**CE**



Velcro fastening



Washable and breathable



4 stays

Ref.	Product	Size	Measurements
902.292	Lumbar back support belt with suspenders	S	86-96 cm x 20,32 cm
902.293	Lumbar back support belt with suspenders	M	96-106 cm x 20,32 cm
902.294	Lumbar back support belt with suspenders	L	106-120 cm x 20,32 cm
902.295	Lumbar back support belt with suspenders	XL	120-135 cm x 20,32 cm
902.275	Strapless lumbar back support belt	S	86-96 cm x 20,32 cm
902.286	Strapless lumbar back support belt	M	96-106 cm x 20,32 cm
902.289	Strapless lumbar back support belt	L	106-120 cm x 20,32 cm
902.290	Strapless lumbar back support belt	XL	120-135 cm x 20,32 cm

## Body protection: ergonomic line

<b>Standard and certification</b>	CE Manufactured according to basic health and safety requirements specified in Regulation (EU) 2016/425																											
<b>Applications</b>	For use in jobs with sedentary postures, sedentary activities, for jobs with air draughts, loading and handling movements, pain due to weather changes, lumbago, rheumatism, and whenever some extra support is needed to ease lower back discomfort.																											
<b>Conservation Storage - Expiry</b>	Store in its packaging in a cool, well-ventilated and dry place. May deteriorate with use. Inspect regularly and replace when very damaged or worn.																											
<b>Directions Use</b>	Wash with tepid water and soap and a cloth or sponge. Do not machine wash, or use aggressive soap. Dry and keep in its box. For personal use. Should not be used by several people, even if they are washed well. Should not be in contact with the skin.																											
<b>Presentation</b>	Individually packed: 1 unit. 50 units per carton.																											
<b>Bar code</b>	<table border="1"> <thead> <tr> <th></th> <th>BAR CODE</th> <th>CARTON BAR CODE</th> </tr> </thead> <tbody> <tr> <td>902.292</td> <td>8423173110338</td> <td>18423173110335</td> </tr> <tr> <td>902.293</td> <td>8423173110345</td> <td>18423173110342</td> </tr> <tr> <td>902.294</td> <td>8423173110352</td> <td>18423173110359</td> </tr> <tr> <td>902.295</td> <td>8423173110369</td> <td>18423173110366</td> </tr> <tr> <td>902.275</td> <td>8423173110161</td> <td>18423173110168</td> </tr> <tr> <td>902.286</td> <td>8423173110277</td> <td>18423173110274</td> </tr> <tr> <td>902.289</td> <td>8423173110307</td> <td>18423173110304</td> </tr> <tr> <td>902.290</td> <td>8423173110314</td> <td>18423173110311</td> </tr> </tbody> </table>		BAR CODE	CARTON BAR CODE	902.292	8423173110338	18423173110335	902.293	8423173110345	18423173110342	902.294	8423173110352	18423173110359	902.295	8423173110369	18423173110366	902.275	8423173110161	18423173110168	902.286	8423173110277	18423173110274	902.289	8423173110307	18423173110304	902.290	8423173110314	18423173110311
	BAR CODE	CARTON BAR CODE																										
902.292	8423173110338	18423173110335																										
902.293	8423173110345	18423173110342																										
902.294	8423173110352	18423173110359																										
902.295	8423173110369	18423173110366																										
902.275	8423173110161	18423173110168																										
902.286	8423173110277	18423173110274																										
902.289	8423173110307	18423173110304																										
902.290	8423173110314	18423173110311																										

